

Innovation and choices  
from your trusted brand

**AnchorFast™**  
Oral Endotracheal Tube Fasteners

Image represents an actor who is not actually intubated or included in clinical study.


**AnchorFast**  
Oral Endotracheal Tube Fastener




**AnchorFast SlimFit**  
Oral Endotracheal Tube Fastener



**AnchorFast Guard**  
Oral Endotracheal Tube Fastener



**AnchorFast Guard Select**  
Oral Endotracheal Tube Fastener



# Our Family of **AnchorFast** Oral Endotracheal Tube Fasteners

Every patient dependent on a ventilator for respiration is vulnerable. Intubated individuals face many risks, and in particular, ventilator-associated events, which can result in longer lengths of stay and increased hospital costs.<sup>1</sup> AnchorFast oral endotracheal tube fasteners combine a number of features that provide critical-care clinicians ease of access to the oral cavity while securing tubes in intubated individuals. AnchorFast oral endotracheal tube fasteners are used in thousands of hospitals around the world; at Hollister, customer feedback drives ongoing innovation, and our product line is continuing to evolve to meet your needs.

## **AnchorFast**

Oral Endotracheal Tube Fastener

**Helps prevent the formation of lip ulcers through the use of the gliding tube shuttle while easing access to the oral cavity**


**Gliding tube shuttle** allows easy side-to-side repositioning of the ET tube for oral care, while the tube remains secure

**Nonabsorbent upper lip stabilizer**

designed to help keep the track from contacting the face

**Adjustable neck strap**

helps provide secondary securement and is replaceable



**ET tube wrap with one-click security clamp**

locks the ET tube in place

**Skin barrier pads**

fasten the device firmly in place. The unique pad formulation is designed for easy application



**Help address intubated patient risks with confidence and offer the right combination of features for your patient's needs with AnchorFast oral endotracheal tube fasteners.**

## AnchorFast SlimFit™

Oral Endotracheal Tube Fastener


Designed to accommodate smaller faces

**Two-piece adjustable neck strap** designed for secondary adjustment at an accessible area on the patient's head

**Gray track** that differentiates from other AnchorFast products

**Thumb pad** with increased surface area for finger or thumb placement when opening

**Nonabsorbent upper lip stabilizer** with increased surface area




## AnchorFast Guard™

Oral Endotracheal Tube Fastener

Designed with an attached bite block

**Integrated tube protection sleeve** helps prevent occlusion without interfering with routine care. Unlike separate bite blocks, it's there *before* it's needed




## AnchorFast Guard Select™

Oral Endotracheal Tube Fastener

Accommodates standard ET tubes and ET tubes with subglottic suctioning

**Modified tube protection sleeve** secures both standard ET tubes from 5.0 to 9.0 mm inner diameter and ET tubes with subglottic suction lumen from 6.0 to 8.0 mm inner diameter



## The AnchorFast oral endotracheal tube fastener product line

	Stock No.	Quantity/ Box	Product Description	Gliding tube shuttle	Nonabsorbent upper-lip stabilizer	ET tube wrap with one-click security clamp	Adjustable neck strap	Skin-barrier pads	Gray track	Two-piece adjustable neck strap	Thumb pad with increased surface area	Nonabsorbent upper-lip stabilizer with increased surface area contact	Integrated tube-protection sleeve	Secures standard and ET tubes with integrated subglottic suction
	9799	12	<b>AnchorFast</b> Oral Endotracheal Tube Fastener Secures ET Tubes ranging from 5.0–10.0 mm	●	●	●	●	●						
	9787	12	<b>AnchorFast SlimFit</b> Oral Endotracheal Tube Fastener Secures ET Tubes ranging from 5.0–8.0 mm	●	●	●	●	●	●	●	●	●		
	9800	12	<b>AnchorFast Guard</b> Oral Endotracheal Tube Fastener Secures ET Tubes ranging from 5.0–10.0 mm	●	●	●	●	●					●	
	9700	12	<b>AnchorFast Guard Select</b> Oral Endotracheal Tube Fastener Secures ET Tubes ranging from 5.0–9.0 mm Secures ET Tubes with integrated subglottic suction lumen ranging from 6.0–8.0 mm	●	●	●	●	●					●	●
	9788	4	<b>AnchorFast</b> Neck Strap Replacement neck band straps for AnchorFast, AnchorFast Guard and AnchorFast Guard Select oral endotracheal tube fasteners											

For more information or to request a sample, visit [hollister.com/anchorfast](http://hollister.com/anchorfast)



**Hollister Incorporated**  
2000 Hollister Drive  
Libertyville, Illinois 60048  
1.800.323.4060

[www.hollister.com](http://www.hollister.com)

<sup>1</sup> Restrepo M, et al. "Economic Burden of Ventilator-Associated Pneumonia Based on Total Resource Utilization." *Infection Control and Hospital Epidemiology*, 2010; 31(5):509-515.

USA Rx Only. Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product. Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

The Hollister logo, AnchorFast, AnchorFast Guard, AnchorFast Guard Select and AnchorFast SlimFit are trademarks of Hollister Incorporated. ©2021 Hollister Incorporated.

923371-0121

